



## SÖHNGEN® ESCORT Skater - Professional Rescue Trolley

The new Escort Skater is an amalgamation of the classic emergency rescue bag with a dynamic trolley system and is designed to easily carry a large assortment of emergency medical products quickly and with the minimum of effort - on wheels or by hand!



ESCORT Skater  
Industry



ESCORT Skater  
Modul Sport



ESCORT Skater  
Modul School XXL



ESCORT Skater  
Modul Senior Citizens



ESCORT Skater  
Modul First Responder



ESCORT Skater  
Emergency Medical Care

Size: 670 x 370 x 350 mm.

Manufactured from durable textured polyvinyl material and finished in a resilient polypropylene coating, the Skater is water withstanding, dirt repellent and easy to clean and maintain after utilisation.

An extensive series of reflection striping on all four sides warrants a high visibility during night time deployment.

Two skate wheels positioned on one end of the robust undercarriage ensure a speedy and safe transport of the skater and its sensitive medical contents. The undercarriage contains two stabilisers to steady the trolley whilst stationary modus .

Opening the bag reveals a large internal storage compartment that is divided into sections using adjustable Velcro based dividers, one long horizontal and two short vertical. These can be adjusted to meet the specific needs and requirements of the user.

Three large module pouches (in red, blue and yellow with transparent covers) are located on the internal lid roof, whilst externally; the lid contains an extra storage compartment that includes two series of rubber lining and loops for the storage of smaller equipment.

The long side compartment, with its series of integrated rubber loops, caters for the storage of diagnostic instruments and other smaller equipment. It opens to a 45° angle.

The head side compartment contains 4 transparent pockets designed for the storage of infusion bags. The lid cover can be removed (with infusion contents) and taken to the accident scene quickly – contains a clothes hanger type hook for hanging.

Available empty or in a series of alternative fillings, each carefully chosen to match requirements of the particular field or discipline.

04 03 009  
04 06 020

ESCORT Skater Trolley • Standard model with DIN 13169 plus Additional Set  
ESCORT Skater Trolley • Empty

04 03 010  
04 03 011  
04 03 012

ESCORT Skater • Sport  
ESCORT Skater • School XXL  
ESCORT Skater • Senior Citizens

04 03 013  
04 03 014  
04 03 015

ESCORT Skater • First Responder  
ESCORT Skater • Emergency Medical Care  
ESCORT Skater • Industry